


EIKONHA

BEYOND THE BEAUTY

INTRODUCTION


Eikonha is next generation cross linked hyaluronic acid dermal filler with lidocaine designed for face , based on Enhanced Interlinked 3D Matrix (EIM) technology which allows easy injectability, natural looking and long-lasting effect.



EIKONHA

EIM Technology

Enhanced interlinked 3D Matrix (EIM) technology is a state of the art manufacturing technique which provides high covalent bond among HA Molecules in a polysaccharide chain.



Benefits of EIM Technology

3-D hydrophilic matrix network

Allows nutrition, oxygen and hormones to go through freely, and protects regular physiological functions of the tissue.

Promising rheological properties

- Increased surface area
- Good Cohesivity of gel


Higher stability against natural degradation by hyaluronidase

GLOW OUT SKIN



EIKONHA GLOW

HA Concentration: 20mg/ml
Lidocaine Concentration: 0.3%
Volumizing Capacity: +++
Volume: 1 ml
Needles & Cannula: 30½ G needle



Injection Depth: Deep Dermis

Indications :


- Perioral lines
- Periorbital lines
- Superficial facial lines and
- Skin depressions
- As a Mesotherapy

POUT OUT LOUD



EIKONHA POUT

HA Concentration: 20mg/ml
Lidocaine Concentration: 0.3%
Volumizing Capacity: ++
Volume: 1 ml
Needles & Cannula: 30½ G needle



Injection Depth: : Mid -Dermis or
Lip Mucosa

Indications :


- Pouting
- Lip fullness
- Vermilion border
- Oral commissures

**GLAM OUT WITH
ELEGANCE**



EIKONHA GLAM

HA Concentration: 20mg/ml
Lidocaine Concentration: 0.3%
Volumizing Capacity: +++
Volume: 1 ml
Needles & Cannula: 27½ G needle



Injection Depth: Deep Dermis

Indications :

- Nasolabial folds
- Cheek's fullness
- Chin augmentation
- Non -Surgical Rhinoplasty

**FILL OUT WITH
CONFIDENCE**



EIKONHA PROFIL


HA Concentration: 20mg/ml

Lidocaine Concentration: 0.3%

Volumizing Capacity: ++++

Volume: 1 ml

Needles & Cannula: 25G Cannula




Injection Depth: Subcutaneous or Supraperiosteal

Indications :

- Volume creation
- Facial contours
- Hands
- Chin and Cheek
- Jawline
- Temple
- Large deep wrinkles and folds

EIKONHA®

Cross-Linked Injectable Hyaluronic Acid



Duration	6m	6m-8m	8m-12m	12m-18m
Volume	1ml	1ml	1ml	1ml
Injection Depth	Intradermal	Mid-dermis Lips mucosa	Deep Dermis	Subcutaneous Supraperiosteal
HA Concentration (mg/mL)+0.3% lidocaine	20mg/ml	20mg/ml	20mg/ml	20mg/ml
Volumizing Capacity	+	++	+++	++++
Needles & Cannula	30G/2 needle in pack	30G/2 needle in pack	27G/2 needle in pack	25G Cannula (1)/ 22G Needle (1)



AK Pharma, Inc., 9100 S Dadeland Blvd Suite 1500-MBP# 150006 Miami, FL 33156

Web: www.eikonha.com | info@eikonha.com